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2.2.3. Recommended procedure

We propose an improved AMP procedure with the ground-level γ -spectrometer as follows:

- 1) Measure the seawater volume (5-100 liters) and put into a tank with appropriate size.
- 2) pH should be adjusted to be 1.6-2.0 by adding concentrated HNO_3 (addition of 40 ml conc. HNO_3 for 20 litre seawater sample makes pH of sample seawater about 1.6).
- 3) Add CsCl of 0.26g to form an insoluble compound and stir at a rate of 25 litre per minute for several minutes.
- 4) Weigh AMP of 4g and pour it into a tank to disperse the AMP with seawater.
- 5) 1 hour stirring at the rate of 25 litre air per minute.
- 6) Settle until the supernate becomes clear. A settling time is usually 6 hours to

overnight, but no longer than 24 hours.

- 7) Take an aliquot of 50 ml supernate to calculate the amount of the residual caesium in the supernate.
- 8) Loosen the AMP/Cs compound from the bottom of the tank and transfer into a 1-2 litre of beaker, if it is necessary do additional step of decantation.
- 9) Collect the AMP/Cs compound onto 5B filter by filtration and wash the compound with 1M HNO₃
- 10) Dry up the AMP/Cs compound for several days in room temperature
- 11) Weigh the AMP/Cs compound and determine weight yield
- 12) Transfer the AMP/Cs compound into a Teflon tube of 4ml volume and subject to γ -ray spectrometry

2.2.4. *Underground γ -spectrometry*

- 1) the same procedure from step 1) to step 12)
- 2) Dissolve AMP/Cs compound by adding alkali solution
- 3) pH should be adjusted to be ca. 8.1 by adding 2M HCl and adjust the volume of solution ca. 70-100ml.
- 4) Perform precipitation of Cs₂Pt(Cl)₄ to add chloroplatinic acid (1g/5ml D.W) at pH = 8.1 and keep in refrigerator during a half-day.
- 5) Collect the Cs₂Pt(Cl)₄ precipitate onto filter by filtration and wash the compound with solution (pH = 8.1)
- 6) Dry up the Cs₂Pt(Cl)₄ precipitate for several days in room temperature
- 7) Weigh the Cs₂Pt(Cl)₄ precipitate and determine weight yield
- 8) Transfer the Cs₂Pt(Cl)₄ precipitate into a Teflon tube of 4ml volume and subject to underground γ -spectrometry

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